

## Pinpoint accuracy for successful treatment

The new **SEIRIN lasers** are eye-catchers not only because of their appearance. What's more, they combine state-of-the-art medical technology, quick and easy operation and ergonomic design. Find out for yourself.



SEIRIN lasers are **wireless, truly "mobile" lasers** that can be used anywhere. A large, bright **colour display** provides ultimate ease of treatment. Another innovation: the display can be rotated at the push of button. This makes SEIRIN lasers ideally suited **for both right-handed and left-handed users and for any treatment position**. In addition, locating acupuncture points is simplified with an **integrated point finder**. As soon as a point is located, laser treatment can be started automatically.


## Convincing laser technology

Whether as **infrared lasers or red light lasers**, all SEIRIN lasers offer **Nogier and Bahr frequencies** in addition to a **Continuous Wave**. The required power can be freely adjusted in 10 mW intervals. In addition to the parameters mentioned above, the SEIRIN laser with 150 mW offers **Reininger frequencies**. Any other required frequencies can be programmed.

Four exchangeable, maintenance-free **high-performance rechargeable batteries** ensure adequate operating times. Like in a mobile phone, batteries are recharged with a practical **fast charger**. No expensive, bulky charging station is needed. Instead, the unit is placed in an elegant acrylic stand.



## Functions of SEIRIN lasers



## Technical specifications

	50 mW infrared	150 mW infrared	40 mW red light
Color display for left-handed and right-handed operation	x	x	x
Indicators for treatment time, laser power, dose, frequency; charge status	x	x	x
Memory function for last saved values	x	x	x
Integrated point finder	x	x	x
Adjustable laser power	10-50 mW	30-150 mW	10-40 mW
Wavelength	785	808	655
Continuous Wave	x	x	x
Nogier frequencies	x	x	x
Bahr frequencies	x	x	x
Reininger frequencies		x	
Freely programmable frequencies		x	
Waveguide, straight	x	x	x
Connection of other waveguides	x	x	x
2 pairs of high-performance rechargeable batteries	x	x	x
Fast charger	x	x	x
Automatic shut-off when not in use	x	x	x
Laser class	3B	3B	3B
Recommended Retail Price	£1350.00	£1922.00	£1350.00
<b>Harmony Medical Introductory Price</b>	<b>£1300.00</b>	<b>£1865.00</b>	<b>£1300.00</b>

*Free UK delivery – prices exclude UK VAT currently at 17.5%*

All SEIRIN products are CE tested and comply with the quality and safety standards of TÜV, ISO 9002 and WHO.

To place your order, call now on 0800 092 8123

Distributor:

### **Harmony Medical**

629 High Road Leytonstone, London E11 4PA

Telephone: +44 (0)20 85187337

Fax: +44 (0)20 8556 5036

E mail: [info@harmonymedical.co.uk](mailto:info@harmonymedical.co.uk)

Website: <http://www.harmonymedical.co.uk>

