

UROstim 2

NEWEST STIMULATOR FOR THE TREATMENT OF INCONTINENCE

The UROstim 2 has been designed to improve the strength of the pelvic floor and sphincter muscle and to treat pain in the gynecological field. The UROstim 2 includes 10 preset programs and 10 free programs for all kinds of incontinence treatment:

- Urge incontinence
- Stress incontinence
- Mixed incontinence
- Anal incontinence

Further the UROstim 2 includes special programs for the treatment of vesical instability and anal incontinence by stimulating the nervous tibialis posterior with self-adhesive electrodes. The parameters in those programs have been validated by several studies showing very good results.

INTEGRATED RECHARGEABLE BATTERY

The UROstim 2 works with an integrated rechargeable battery which reduces the long time cost for both patients and therapists.


SUPPLIED WITH THE FOLLOWING ACCESSORIES

- 106351 Cable type 5.15 (2 pieces)
- 283400 Self-adhesive electrodes 50 x 50 mm (4 pieces)
- 602000 9 V battery
- Box for easy transportation
- Operating manual

SEE ALL ACCESSORIES INCLUDING VAGINAL AND ANAL PROBES ON PAGES 22 TO 27



FEATURES



TYPE OF PROGRAM

- | | | |
|--|--------------------------------------|--|
| <input checked="" type="checkbox"/> PAIN | <input type="checkbox"/> MUSCLE STIM | <input type="checkbox"/> DYNAMIC |
| <input type="checkbox"/> HAN | <input type="checkbox"/> DENTAL | <input checked="" type="checkbox"/> UROLOGICAL |
| <input type="checkbox"/> SPORT | <input type="checkbox"/> BEAUTY | <input type="checkbox"/> WELLNESS |
| <input type="checkbox"/> REHABILITATION | <input type="checkbox"/> RELAXATION | <input type="checkbox"/> SPECIAL |

TECHNICAL DATA

Item no.	101453
Output current	100 mA (with 1 kΩ real)
Frequency	0.5 Hz to 120 Hz
Pulse width	60 µs to 300 µs
Power supply	Integrated rechargeable battery