

## SPORÉCUP XTR2

### MUSCLE STIMULATOR FOR SPORTS AND FITNESS

The XTR2 is a TENS/EMS combination unit specially designed for athletes. It includes 30 programs divided in the following categories:

- Recovery
- Relaxation
- Pain therapy
- Sport
- Rehabilitation
- Capillarisation
- Fitness
- Aesthetics

The XTR2 is the newest muscle stimulator in the schwa-medico range. It combines high quality stimulation with modern and friendly design.

### INTEGRATED RECHARGEABLE BATTERY

The XTR2 works with an integrated rechargeable battery which reduces the long time cost for batteries for both patients and therapists.


### SUPPLIED WITH THE FOLLOWING ACCESSORIES

- |        |  |
|--------|--|
| 106351 | Cable type 5.15 (2 pieces)                     |
| 283400 | Self-adhesive electrodes 50 x 50 mm (4 pieces) |
| 101062 | Charger WK 112 for XTR2                        |
|        | Transportation case                            |
|        | Operating manual                               |

### SEE OPTIONAL ACCESSORIES ON PAGES 22 TO 25



### FEATURES



### TYPE OF PROGRAM

- |   |   |  |
|---|---|--|
| <input checked="" type="checkbox"/> PAIN  | <input checked="" type="checkbox"/> MUSCLE STIM | <input type="checkbox"/> DYNAMIC             |
| <input type="checkbox"/> HAN              | <input type="checkbox"/> DENTAL                 | <input type="checkbox"/> UROLOGICAL          |
| <input checked="" type="checkbox"/> SPORT | <input checked="" type="checkbox"/> BEAUTY      | <input checked="" type="checkbox"/> WELLNESS |
| <input type="checkbox"/> REHABILITATION   | <input checked="" type="checkbox"/> RELAXATION  | <input type="checkbox"/> SPECIAL             |

### TECHNICAL DATA

Item no.	104064
Output current	100 mA (with 1 kΩ real)
Frequency	0.5 – 120 Hz
Pulse width	50 – 400 μs
Power supply	Integrated rechargeable battery