

Vaginal and anal probes

To be used with UROstim, URO stim 2, EMP 2 PRO and EasyTrain MYO for the treatment of incontinence (type 5 cables).

VAGINAL PROBE VATT

For use in cases of vaginal convergence and shortening

103042 Length 55 mm, Diameter 20 mm



VAGINAL PROBE V2B

For use in cases of vaginal shortening

103025 Length 58 mm, Diameter 25 mm



VAGINAL PROBE TAM

For use in cases of vaginal shortening

103043 Length 58 mm, Diameter 25 mm



ANAL PROBE TL1R

With stop at the back end for fixation. The thickening at the front prevents the electrode from being ejected.

103044 Length 65 mm, Diameter 12 - 18 mm



ANAL PROBE 12C

The flexible stop at the back end with a positioning bar and the conical end at the front provide permanent attachment during treatment.

103045 Length 64 - 108 mm, Diameter 12 - 28 mm

